

Hi Mel!

Thank you again so so so so much for doing this, you are a gem. Did I mention you're welcome to my Anchor steam? You're welcome to my Anchor Steam, in the fridge and next to the fridge.

Here are dates w/ how I'd probably be doing them; please feel free to modify as you like/see fit, just let me know what you end up doing. It also turns out upon further examination that I fucked up the written lab notebook calendar, so...ignore it? Ish? This is way easier to read.

Relevant days in gray.

P19 day	Date	action
0	Sep 25	Plate new cells,
1	Sep 26 Friday	Switch 5 of 6 of the P19 dishes to 1% MEMalpha with 1 mM RA; the last dish stays on 5% to act as a stock
3	Sep 28 Sunday	Change media on all dishes Trizol lysing of one set of diff cells
6	Oct 1 Wednesday	Change media lyse one dish split stock P19s (maybe) note pass #, date, media type, cell type,
9	Oct 4 Saturday	Change media lyse one dish
	I get back October 5 th ! Yay!	
12	Oct 7 Tuesday	Change media lyse one dish split stock P19s

protocols for Baby Mel!

differentiating P19s

Why?

P19 stem cells differentiate into neuron-like cells when exposed to 1 mM trans-retinoic acid (RA). While RA technically only needs to be added for the first 3 days or so to set cell fate, we keep it on RA the whole time and on low levels of FBS to promote differentiation.

protocol

reagents

- retinoic acid (-20 C freezer, in the corner with trypsin; kept in tinfoil to prevent light exposure and degradation)
- MEMalpha with 1% FBS
- sparse-ish P19, preferably at a low pass number

do

1. Take off any current media (typically 5%) and wash the cells with PBS gently (if this is Day 0 of differentiation, otherwise no washing!). P19s throw off a lot of stuff.
2. Add 1.5 mls MEMalpha with 1% FBS
 - the less FBS, the better they differentiate, but also the faster they die. It's a balance.
3. Add 1.5 uls of the RA to the center of the dish; swirl!
 - a cool thing happens where the RA must be very hydrophobic, so a tiny clear spot will appear where you pipette it in
4. Do this again 3 days later!
 - P19s are delicate infants; please don't wash them if they're in the process of throwing out processes.

RNA extraction from P19s

Why?

Extracting, isolating, and quantiating RNA expression allows us to build a timeline of the "usual" P19 differentiation expression levels, and then compare it to what happens with the introduction of different transcription factors. Step 1: Kill cells, collect cell lysate to later do an RNA extraction en masse.

Do

1. Wash cells in 35mm dish 2x with PBS (wash= add ~1ml PBS and aspirate).
 - unless you're using delicate cells that don't need washing, like delicate differentiating P19s
2. Add 1 mL of TRIZOL Reagent to the dish and pipette up and down to make sure all the cells are really dead and really broken up.
3. Scrape TRIZOL/lysed cells into a 1.5 mL microcentrifuge; scrape! Get ALL of the stuff.
4. Incubate the homogenized samples for 5 minutes at room temperature
 - permits complete dissociation of nucleoprotein complexes.

Store the suspension @ -80 [in the Leeper box, green tape, top shelf next to the other mol. bio stuff] , label with date & any other pertinent information

Splitting P19s

why?

sometimes cells are too confluent and we need less of them! P19s grow fast, so even with a 1:9 split (no counting @ Hampshire), they need to be split every 5-6 days, and fed every three days.

do

reagents

- MEMalpha, with some percent FBS depending on your task. 5% for normal maintenance.
- trypsin
- 1. Aspirate media, rinse with PBS
- 2. Add .5 ml trypsin and wait for cells to shake loose
- if you're not me, you'll actually let them hang out in the incubator and dissociate properly
- 3. Add 1.5 ml media to deactivate trypsin
- 4. add .2 mls cells and 1.8 mls media to each plate for a ~1:9 split